

Sending Data from Microcontroller to Cloud via WiFi (Mobile Hotspot)

1. Introduction

In this project, the **ESP32 microcontroller** reads sensor data (heart rate, SpO₂, temperature) from the patient and sends it to **Firestore Realtime Database**. The data is transmitted using **WiFi**, which can be provided via a **mobile hotspot** if no other internet connection is available.

2. Why Use Mobile Hotspot

- Mobile hotspot provides **temporary internet access** for ESP32.
 - Useful in remote areas or during testing where a router or fixed WiFi is unavailable.
 - ESP32 can connect to the hotspot just like a normal WiFi network.
-

3. ESP32 Connection Process

1. WiFi Credentials:

- ESP32 stores the **SSID (network name)** and **password** of the mobile hotspot.
- When powered on, it searches for the hotspot and connects.

2. Connection Check:

- ESP32 continuously checks whether it is **successfully connected** to the hotspot.
- If connection fails, it can **retry automatically** until connected.

3. Internet Access:

- Once connected, ESP32 gains access to the internet.
 - It can now communicate with **Firestore servers** through the cloud.
-

4. Sending Data to Firestore

- After WiFi is connected, ESP32 uses the **Firestore ESP Client library** to authenticate with Firestore.
 - Data from sensors is read and prepared in the **JSON format**.
 - ESP32 sends the data to specific paths in Firestore (e.g., `patients/P102/heart_rate`).
 - Firestore updates instantly, allowing **real-time monitoring**.
-

5. Data Flow Summary

1. ESP32 reads sensor values.
2. ESP32 connects to mobile hotspot WiFi.

3. WiFi provides internet access to ESP32.
 4. ESP32 authenticates with Firebase using **API Key** and **Database URL**.
 5. Sensor data is sent to Firebase in **JSON tree format**.
 6. Firebase updates **live**, pushing data to connected applications (Android/web dashboard).
-

6. Advantages of Using Mobile Hotspot

- Provides **flexible connectivity** without relying on home or office WiFi.
 - Useful for **testing and field deployment**.
 - Works anywhere a mobile network is available.
 - Allows **real-time data upload** for patient monitoring.
-

7. Considerations

- Mobile hotspot must have a **stable internet connection**.
 - Power consumption on ESP32 may increase if WiFi connection is weak.
 - For multiple patients or long-term monitoring, a **dedicated WiFi network** is more reliable.
-

8. Conclusion

Using a **mobile hotspot** enables the ESP32 to send real-time patient data to Firebase efficiently, making the system flexible and deployable in areas without fixed internet. This setup ensures **live monitoring, instant updates, and secure cloud storage** for patient health information.